

Straight Sharp Cannulae



Straight Blunt Cannulae



Curved Sharp Cannulae




ESENTEC™

Conventional RFA

A CONVENTIONAL LESION OPTION

The ESENTEC™ line offers RFA cannulae and probes in a variety of gauge sizes and lengths. Available in a variety of tip styles such as straight sharp, curved sharp, blunt, and curved blunt. ESENTEC offers options with a radiopaque band to aid in fluoroscopic visualization.



IN CONVENTIONAL RF THE LESIONS ARE ELIPTICAL†

With Conventional RFA, the lesion projects distally a minimal distance beyond the tip, resulting in an elliptical shape.†

AVANOS*



Call 1-844-4AVANOS
(1-844-428-2667)
in the United States and Canada

Visit avanospainmanagement.com for the product catalog and ordering information.

OUR MISSION

In the face of the global opioid epidemic, the AVANOS Pain Management & Recovery portfolio provides non-opioid solutions tailored for the management and treatment of chronic and post-surgical pain.

AVANOS continues to push toward a growing portfolio of opioid-sparing therapies. This includes HA injections and radiofrequency nerve ablation for chronic pain, cryocompression therapy for injuries, and post-surgical pain pumps for recovery. Each product is designed to **help patients get back to the things that matter.**[†]



[†] Claim data is on file at Avanos