


AVANOS*

NOPAIN ACT PRODUCT CATALOG



ON-Q*
PAIN RELIEF SYSTEM

ambiT
AMBULATORY INFUSION THERAPY

ON-Q* ELASTOMERIC PUMPS

PRODUCT CODE	FILL VOLUME X FLOW RATE	MAX FILL VOLUME	BOLUS	PKG/ CASE
ON-Q* PUMP WITH SELECT-A-FLOW*				
SAF01	270 ml, 1-7 ml/hr	335 ml	–	5
CB004	400 ml, 2-14 ml/hr	550 ml	–	5
CB6004	600 ml, 2-14 ml/hr	750 ml	–	5
CB4007	Dual, 400 ml, 1-7 ml/hr per site	550 ml	–	5
CB6007	Dual, 600 ml, 1-7 ml/hr per site	750 ml	–	5
ON-Q* PUMP WITH SELECT-A-FLOW* & ONDEMAND*				
CB005	400 ml, 1-7 ml/hr with 5 ml bolus (60-minute refill)	550 ml	5 ml/60 minutes	5
CB006	400 ml, 2-14 ml/hr with 5 ml bolus (30-minute refill)	550 ml	5 ml/30 minutes	5
ON-Q* PUMP WITH ONDEMAND*				
CB003	400 ml, 5ml/hr with 5 ml bolus (60-minute refill)	550 ml	5 ml/60 minutes	
ON-Q* FIXED FLOW PUMP - PUMP ONLY MODELS				
P100X2	100 ml x 2 ml/hr	125 ml	–	5
P100X2D	Dual, 100 ml x 1 ml/hr per site	125 ml	–	5
P270X2	270 ml x 2 ml/hr	335 ml	–	5
P270X4	270 ml x 4 ml/hr	335 ml	–	5
P270X5	270 ml x 5 ml/hr	335 ml	–	5
P270X4D	Dual, 270 ml x 2 ml per site	335 ml	–	5
P400X4	400 ml x 4 ml/hr	550 ml	–	5
P400X8	400 ml x 8 ml/hr	550 ml	–	5
P400X4D	Dual, 400 ml x 2 ml per site	550 ml	–	5
CB001	400 ml x 5 ml/hr	550 ml	–	5
CB002	400 ml x 10 ml/hr	550 ml	–	5

PRODUCT CODE	FILL VOLUME X FLOW RATE	MAX FILL VOLUME	PKG/ CASE
PUMP WITH ON-Q* 1 IN. SILVERSOAKER* CATHETER			
PM032-A	100 ml x 2 ml/hr	125 ml	5
PM036-A	100 ml x 1 ml/hr	125 ml	5
PUMP WITH ON-Q* 2.5 IN. SILVERSOAKER* CATHETER			
PM012-A	100 ml x 2 ml/hr	125 ml	5
PM013-A	270 ml x 2 ml/hr	335 ml	5
PM014-A	270 ml x 5 ml/hr	335 ml	5
PM015-A	Dual, 270 ml x 2 ml/hr per site	335 ml	5
PM018-A	Dual, 400 ml x 2 ml/hr per site	335 ml	5
PUMP WITH ON-Q* 5 IN. SILVERSOAKER* CATHETER			
PM022-A	100 ml x 2 ml/hr	125 ml	5
PM023-A	270 ml x 2 ml/hr	335 ml	5
PM024-A	270 ml x 5 ml/hr	335 ml	5
PM025-A	Dual, 270 ml x 2 ml/hr per site	335 ml	5
PM026-A	400 ml x 5 ml/hr	550 ml	5
PM028-A	Dual, 400 ml x 2 ml/hr per site	550 ml	5
PUMP WITH ON-Q* 10 IN. SILVERSOAKER* CATHETER			
PM045-A	Dual, 270 ml x 2 ml/hr per site	550 ml	5
PM046-A	400 ml x 5 ml/hr	550 ml	5
PM048-A	Dual, 400 ml x 2 ml/hr per site	550 ml	5

ambIT® ELECTRONIC PAIN CONTROL SYSTEM PRODUCTS

PRODUCT CODE	DESCRIPTION	UNITS PER SHIPPING BOX
ambIT*+ PROGRAM SOLUTIONS		
220575	ambIT + Reusable w/Filter Spike Cassette - Loaned PIB-IL 220527 Pump Spike Cassette with Filter, Carrying Pouch, Two (2) "AA" Batteries, Button Guard w/o access, and Patient Manual	5
220578	ambIT + Disposable w/Filter Spike Cassette - PIB-PCA Pump Spike Cassette with Filter, Carrying Pouch, Two (2) "AA" Batteries, Button Guard w/o access, and Patient Manual	5
ambIT* REUSABLE SOLUTIONS		
STANDARD KIT (CASSETTE + PUMP ACCESSORIES)		
220316	ambIT Pain Control Kit - Basic with Filterless Spike Cassette Cassette Filterless Spike, Carrying Pouch, Two (2) "AA" Batteries, Setup Instructions, and Patient Manual	5
220463	ambIT Pain Control Kit - Basic with Filter Spike Cassette Cassette Filter Spike, Button Guard, Carrying Pouch, Two (2) "AA" Batteries, Setup Instructions, and Patient Manual	5
220466	ambIT Pain Control Kit - Basic with Male Luer Filter Cassette and Shipper Cap Cassette Filter Male Luer, Shipper Cap, Carrying Pouch, Two (2) "AA" Batteries, Setup Instructions, and Patient Manual	5
220537	ambIT Pain Control Kit - Basic with Male Luer Filter Cassette and Button Guard Cassette Filter Male Luer, Two AA Batteries, Carrying Pouch, Pump Insert, Patient Manual, and Button Guard w/o Button Access	5
REUSABLE PUMP		
220244	ambIT Continuous Pump Uses: Chemotherapy, Anesthesia, Hydration, and Other Continuous Rate Therapies; Pump Color: Bone Gray	1
220262	ambIT PCA Pump Uses: Post-Operative Local Pain Management, Regional Nerve Blocks, PCA (patient controlled analgesia) Narcotics; Pump Color: Bone Gray; Includes Port for Optional Remote Bolus Button	1
220437	ambIT Military PCA Pump Uses: Approved for all military pain missions, including IV PCA, Epidural, and Continuous Nerve Infusion; Pump Color: International Orange	1
220501	ambIT PIB-PCA Pump Uses: Post-Operative Local Pain Management, Regional Nerve Blocks, PCA (patient controlled analgesia), PIB (Patient Intermittent Bolus), Epidurals, Narcotics Pump Color: Aqua; Includes Port for Optional Remote Bolus Switch	1
220527	ambIT PIB IL Pump Uses: Post-Operative Local Pain Management, Regional Nerve Blocks, PCA (patient controlled analgesia), PIB (Patient Intermittent Bolus), Epidurals, Narcotics Pump Color: Aqua; Includes Port for Optional Remote Bolus Switch	1

PRODUCT CODE	DESCRIPTION	UNITS PER SHIPPING BOX
ambIT DISPOSABLE SYSTEM SOLUTIONS (DISPOSABLE PUMP + CASSETTE + PUMP ACCESSORIES)		
220342	ambIT Preset Pain Control System - Basic with Filterless Male Luer Cassette ambIT PreSet Pump, Cassette Filterless Male Luer, Carrying Pouch, Two (2) "AA" Batteries, Setup Instructions, and Patient and Clinician Manuals	5
220343	ambIT Preset Pain Control System - Basic with Filterless Spike Cassette ambIT PreSet™ Pump, Cassette Filterless Spike, Carrying Pouch, Two (2) "AA" Batteries, Setup Instructions, and Patient/Clinician Manuals	5
220462	ambIT Preset Pain Control System - Basic with Filter Spike Cassette ambIT PreSet Pump, Cassette Filter Spike, Carrying Pouch, Two (2) "AA" Batteries, Setup Instructions, and Patient/Clinician Manuals	5
220563	ambIT Preset Pain Control System - Filter Spike Cassette and Button Guard ambIT PreSet Pump, Spike Cassette with Filter, Carrying Pouch, Two (2) "AA" Batteries, Button Guard, Setup Instructions, and Patient/Clinician Manuals	5
220573	ambIT PIB/PCA v1.5 Pain Control System - Luer Cassette ambIT PIB/PCA v1.5 Pump, Male Luer Cassette with Filter, Carrying Pouch, Two (2) "AA" Batteries, Button Guard, Setup Instructions, and Patient/Clinician	5
220574	ambIT PIB/PCA v1.5 Pain Control System - Spike Cassette ambIT PIB/PCA v1.5 Pump, Spike Filter Cassette, Carrying Pouch, Two (2) "AA" Batteries, Button Guard, Setup	5
ambIT ACCESSORIES		
220320	Replacement Battery Cap	1
220409	ambIT Pump Carrying Pouch	1
220444	Protective Shipping Cap	1

COPY-08401 Rev 2

INCISIONAL ACCESSORIES

SOAKER* PRODUCT CODE	SILVERSOAKER* PRODUCT CODE	INFUSION LENGTH	PKG/CASE
INCISIONAL CATHETERS			
PM010	PM010-A	Soaker* 2.5 in (6.5 cm), T-Peel Needle 3.25 in	5
PM020	PM020-A	Soaker* 5 in (12.5 cm), T-Peel Needle 6 in	5
PM030	PM030-A	Soaker* 1 in (2.5 cm), T-Peel Needle 3.25 in	5
PM040	PM040-A	Soaker* 10 in (25 cm), T-Peel Needle 8 in	5
-	PM050-A	Soaker* 7.5 in (19.1 cm), T-Peel Needle 8 in	5
PRODUCT CODE	DESCRIPTION		PKG/CASE
INTRODUCER NEEDLES			
ACC01	T-Peel Needle 3.25 in		10
ACC02	T-Peel Needle 6 in		10
ACC03	T-Peel Needle 8 in		10

PRODUCT CODE	DESCRIPTION	PKG/CASE
TUNNELING SYSTEM		
ACC07S	Tunneler Sheath 11 GA x 8 in	10
ACC08S	Tunneler Sheath 17 GA x 8 in	10
ACC09S	Tunneler Sheath 11 GA x 12 in	10
ACC10S	Tunneler Sheath 16 GA x 12 in	10
T17X5	Disposable Tunneler and Sheaths 17 GA x 5 in	5
T17X8	Disposable Tunneler and Sheaths 17 GA x 8 in	5
T17X10	Disposable Tunneler and Sheaths 17 GA x 10 in	5
T16X12	Disposable Tunneler and Sheaths 16 GA x 12 in	5
T11X10	Disposable Tunneler and Sheaths 11 GA x 10 in	5
T11X8	Disposable Tunneler and Sheaths 11 GA x 8 in	5
T11X12	Disposable Tunneler and Sheaths 11 GA x 12 in	5
FT11X12	Disposable Fluid Tunneler and Sheaths 11 GA x 12 in	5
FT11X8	Disposable Fluid Tunneler and Sheaths 11 GA x 8 in	5
FT16X12	Disposable Fluid Tunneler and Sheaths 16 GA x 12 in	5
FT16X8	Disposable Fluid Tunneler and Sheaths 16 GA x 8 in	5
PRODUCT CODE	DESCRIPTION	PKG/CASE
CATHETER PLACEMENT ACCESSORY		
MI22	MACB™ Introducer - Dual Slotted Cannula & Cannulated Trocha	10

COPY-08401 Rev 2

NEEDLES & CATHETER SETS

PRODUCT CODE	NEEDLE						CATHETER		PKG/ CASE
	TIP: B F FACET BEVEL T TUOHY	SL SLEEVE INSULATED SP SPRAY INSULATED E ECHOGENIC	GAUGE	USABLE LENGTH (MM/IN)	STIMULATING CABLE	INFUSION SET	GAUGE	LENGTH (IN)	
ON-Q* QUIKBLOC* OVER-THE-NEEDLE CATHETER SETS									
QB10016SG	F	E	20	100/4		•	16	3	5
QB10016SGC	F	E	5	100/4	•	•	16	3	5
QB15016SG	F	E	5	150/6		•	16	5	5
QB15016SGC	F	E	5	150/6	•	•	16	5	5
PROLONG* CONTINUOUS NEEDLE/CATHETER SETS									
PL18050TGC	T	SP	18	50/2	•	•	21	24	10
PL18100TGC	T	SP	18	100/4	•	•	21	24	10
PL18150TGC	T	SP	18	150/6	•	•	21	24	10
PL19050GCA	B	SP	19	50/2	•	•	21	24	10
PL19100GC	B	SP	19	100/4	•	•	21	24	10
PL19150GC	B	SP	19	150/6	•	•	21	24	10

PRODUCT CODE	NEEDLE						CATHETER		PKG/ CASE
	TIP: B F FACET BEVEL T TUOHY	SL SLEEVE INSULATED SP SPRAY INSULATED E ECHOGENIC	GAUGE	USABLE LENGTH (MM/IN)	STIMULATING CABLE	INFUSION SET	GAUGE	LENGTH (IN)	
ECHOLONG* ECHOGENIC CONTINUOUS NEEDLE/CATHETER SETS									
ELL18100SGC	F	E	18	100/4	•	•	21	24	10
ELL18050TG	T	E	18	50/2		•	21	24	10
ELL18100SGC	F	E	18	100/4	•	•	21	24	10
ELL18050TG	T	E	18	50/2		•	21	24	10
ELL18050TGC	T	E	18	50/2	•	•	21	24	10
ELL18100TG	T	E	18	100/4		•	21	24	10
ELL18100TGC	T	E	18	100/4	•	•	21	24	10

T-BLOC* CONTINUOUS TRAYS

CONTINUOUS TRAYS

COMPONENTS

	TBT02050ST 2" STIM/ECHO NEEDLE CONTINUOUS TRAY WITH DRUGS	TBT02100ST 4" STIM/ECHO NEEDLE CONTINUOUS TRAY WITH DRUGS	TBT02150ST 6" STIM/ECHO NEEDLE CONTINUOUS TRAY WITH DRUGS	TBT03050ST 2" STIM/ECHO NEEDLE CONTINUOUS TRAY WITHOUT DRUGS	TBT03100ST 4" STIM/ECHO NEEDLE CONTINUOUS TRAY WITHOUT DRUGS	TBT02050T 2" ECHO NEEDLE CONTINUOUS TRAY WITH DRUGS	TBT02100T 4" ECHO NEEDLE CONTINUOUS TRAY WITH DRUGS	TBT02150T 6" ECHO NEEDLE CONTINUOUS TRAY WITH DRUGS	TBT03100T 4" ECHO NEEDLE CONTINUOUS TRAY WITHOUT DRUGS
18 GA x 2" (50 mm) EchoLong* Echogenic Tuohy Needle, Stimulating	•			•					
18 GA x 4" (100 mm) EchoLong* Echogenic Tuohy Needle, Stimulating		•			•				
18 GA x 6" (150 mm) EchoLong* Echogenic Tuohy Needle, Stimulating			•						
18 GA x 2" (50 mm) EchoLong* Echogenic Tuohy Needle, Non-Stimulating						•			
18 GA x 4" (100 mm) EchoLong* Echogenic Tuohy Needle, Non-Stimulating							•		•
18 GA x 6" (150 mm) EchoLong* Echogenic Tuohy Needle, Non-Stimulating								•	
20 GA x 24" SPIROL® Catheter, Closed Tip, 3 ports	•	•	•	•	•	•	•	•	•
Catheter Connector	•	•	•	•	•	•	•	•	•
Catheter Securement Device, LockIt Plus®	•	•	•	•	•	•	•	•	•
Connector Securement, Grip-Lok®	•	•	•	•	•	•	•	•	•
ChlorPrep® 3 ml, Orange Tint, 1 or 2	2	2	2	2	2	2	2	2	2
UltraSound Gel	•	•	•	•	•	•	•	•	•
Ultrasound Probe Cover, Elastic	•	•	•	•	•	•	•	•	•
Syringe, 5 cc with Luer Lock	•	•	•	•	•	•	•	•	•
Syringe, 20 cc with Luer Lock	•	•	•	•	•	•	•	•	•
Local Injection needle 18 GA x 1.5" (38 mm)	•	•	•	•	•	•	•	•	•
Local Injection needle 25 GA x 1.5" (38 mm)	•	•	•	•	•	•	•	•	•
Filling Extension Set, 24"	•	•	•	•	•	•	•	•	•
Filling Extension Set, 6"	•	•	•	•	•	•	•	•	•
Lidocaine, 1% Concentration 5 ml Ampule	•	•	•	•	•	•	•	•	•
CSR Wrap for Sterile Field (30" x 30")	•	•	•	•	•	•	•	•	•
Ultrasound Drape (24" x 32" with Perforation and 5.5" x 5.5" Fenestration) with Adhesive Around Aperture	•	•	•	•	•	•	•	•	•
Needle Stick/Disposal Pad	•	•	•	•	•	•	•	•	•
Tegaderm® Dressings 10 cm x 12 cm, 2	•	•	•	•	•	•	•	•	•
Steri-Strip®, 6	•	•	•	•	•	•	•	•	•
4" x 4" Gauze Pads, 8	•	•	•	•	•	•	•	•	•
Filter Straw, 5 Micron	•	•	•	•	•	•	•	•	•
Benzoin, ⅓ cc Vial	•	•	•	•	•	•	•	•	•
MRI Label	•	•	•	•	•	•	•	•	•
Medication ID Labels	•	•	•	•	•	•	•	•	•

T-BLOC* CATHETER SETS & PREP TRAYS

NEEDLE/CATHETER SETS

PREP TRAYS

COMPONENTS

	TB050ST 2" STIM/ECHO NEEDLE NEEDLE/ CATHETER SET	TB100ST 4" STIM/ECHO NEEDLE NEEDLE/ CATHETER SET	TB150ST 6" STIM/ECHO NEEDLE NEEDLE/ CATHETER SET	TB050T 2" ECHO NEEDLE NEEDLE/ CATHETER SET	TB100T 4" ECHO NEEDLE NEEDLE/ CATHETER SET	TB150T 6" ECHO NEEDLE NEEDLE/ CATHETER SET	TBT02 PREP TRAY WITH DRUGS	TBT03 PREP TRAY WITHOUT DRUGS
18 GA x 2" (50 mm) EchoLong* Echogenic Tuohy Needle, Stimulating	•							
18 GA x 4" (100 mm) EchoLong* Echogenic Tuohy Needle, Stimulating		•						
18 GA x 6" (150 mm) EchoLong* Echogenic Tuohy Needle, Stimulating			•					
18 GA x 2" (50 mm) EchoLong* Echogenic Tuohy Needle, Non-Stimulating				•				
18 GA x 4" (100 mm) EchoLong* Echogenic Tuohy Needle, Non-Stimulating					•			
18 GA x 6" (150 mm) EchoLong* Echogenic Tuohy Needle, Non-Stimulating						•		
20 GA x 24" SPIROL® Catheter, Closed Tip, 3 ports	•	•	•	•	•	•		
Catheter Connector	•	•	•	•	•	•		
Catheter Securement Device, LockIt Plus®	•	•	•	•	•	•		
Connector Securement, Grip-Lok®	•	•	•	•	•	•		
ChlorPrep® 3 ml, Orange Tint, 1 or 2							2	2
UltraSound Gel							•	•
Ultrasound Probe Cover, Elastic							•	•
Syringe, 5 cc with Luer Lock							•	
Syringe, 20 cc with Luer Lock							•	•
Local Injection needle 18 GA x 1.5" (38 mm)							•	•
Local Injection needle 25 GA x 1.5" (38 mm)							•	•
Filling Extension Set, 24"							•	•
Filling Extension Set, 6"							•	
Lidocaine, 1% Concentration 5 ml Ampule							•	
CSR Wrap for Sterile Field (30" x 30")							•	•
Ultrasound Drape (24" x 32" with Perforation and 5.5" x 5.5" Fenestration) with Adhesive Around Aperture							•	•
Needle Stick/Disposal Pad							•	•
Tegaderm® Dressings 10 cm x 12 cm, 2		•	•	•	•	•	•	•
Steri-Strip®, 6							•	•
4" x 4" Gauze Pads, 8							•	•
Filter Straw, 5 Micron							•	
Benzoin, ⅓ cc Vial							•	•
MRI Label		•	•	•	•	•		
Medication ID Labels							•	•

ADDITIONAL ACCESSORIES COVERED UNDER THE NOPAIN ACT

Codes include the ambIT or ON-Q pump and all disposable supplies required for catheter insertion

Catheter
Prep tray
Ultrasound sleeve and gel
Surgical drape
Skin adhesive
Dressings

For more information on the NOPAIN Act, please click here: [avano.com/NOPAIN](https://www.avanos.com/NOPAIN), and/or contact our Market Access and Reimbursement Team at reimbursement@avano.com.

1-844-4AVANOS



There are inherent risks in all medical devices. Please refer to the product labeling for Indications, Cautions, Warnings and Contraindications. Failure to follow the product labeling could directly impact patient safety. Physician is responsible for prescribing and administering medications per instructions provided by the drug manufacturer. Refer to www.avanos.com for additional product safety Technical Bulletins.

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
© 2018-2026 AVNS. All rights reserved. COPY-08401 Rev 2