




# NG/NI TUBE RETENTION SYSTEM OPTIONS

**CORGRIP\* | CORGRIP\*SR**




**AVANOS\***

# CORGRIP\* NG/NI Feeding Tube Retention System



PRINTED CENTIMETER MARKINGS



RETRIEVAL CATHETERS

BIOCOMPATIBLE, POLYESTER UMBILICAL TAPE



LUBRICANT




SLOT-SIDE-LOCK COLOR-CODED CLIPS



ORANGE OPENING DEVICE


Product Code	Clamp Size (FR)	QTY / Case
25-008	8	10
25-010	10	10
25-012	12	10
25-014	14	10
25-016	16	10
25-018	18	10

# CORGRIP® SR Enhances Visibility and Flexibility



## New Feature

**RETRIEVAL CATHETER** has a hollow tip that allows for flexibility around the vomer bone†



**PRINTED CENTIMETER MARKINGS** on the catheters help facilitate proper insertion depth†

## New Feature

**CORGRIP® SR**, featuring a **MONOFILAMENT THREAD**, is the nasal bridle system offering numerous advantages vs. tape<sup>1-5</sup>


## New Feature

**YELLOW TIP** designed for enhanced visibility of the magnet†

**LUBRICANT**



**SLOT-SIDE-LOCK COLOR-CODED CLIPS** help facilitate French size differentiation†



**ORANGE OPENING DEVICE** designed for use with CORGRIP® SR†



Product Code	Clamp Size (FR)	QTY / Case
26-008	8	10
26-010	10	10
26-012	12	10
26-014	14	10
26-016	16	10
26-018	18	10

# AVANOS\*

Enteral Feeding

Call 1-844-4AVANOS  
(1-844-428-2667)  
in the United States and Canada

For more information about all AVANOS\* products and resources, visit our websites:

For Patients: **tubefed.com**  
BY AVANOS

For Clinicians: [www.avanosmedicaldevices.com](http://www.avanosmedicaldevices.com)

**At AVANOS, we share your commitment to helping enteral feeding patients get the nutrition they need.** This drives us to continuously improve and grow our offering to provide your patients with products and services to get them back to the things that matter.

A broad enteral feeding portfolio from neonates and pediatrics to adults, along with the dedicated education, service, and support you and your patients need, is what makes **AVANOS enteral feeding a vital lifeline.**<sup>†</sup>



There are inherent risks in all medical devices. Please refer to the product labeling for indications, cautions, warnings and contraindications by visiting <https://eifu.avanos.com/AVA/all>. Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region. Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.

<sup>†</sup> Claim data is on file at Avanos. 1 Lynch A, et al. A systematic review of the effectiveness and complications of using nasal bridles to secure nasoenteral feeding tubes. Aust J Otolaryngol. 2018;1(8):1-11. 2 Gunn SR, et al. Use of a Nasal Bridle Prevents Accidental Nasoenteral Feeding Tube Removal. JPEN Parenter Enteral Nutr. 2009;33(1):50-54. 3 Curram KC. Nasal Bridle: Married to Your Tube. Pract Gastroenterol. January 2011;91:27-34. 4 Allan K, et al. Observation of inadvertent tube loss in ICU: effect of nasal bridles. Br J Nurs. 2019 Oct 10;28(18):1170-1174. 5 Burns SM, et al. Comparison of nasogastric tube securing methods and tube types in medical intensive care patients. Am J Crit Care. 1995 May;4(3):198-203.

Actor portrayals - not an actual patient. | Caution: Federal (U.S.) law restricts this device to sale by, or on the order of a Physician. | Rx Only

[www.avanosmedicaldevices.com](http://www.avanosmedicaldevices.com)

\*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates. © 2018-2025 AVNS. All rights reserved. COPY-08258